
ATTENTION:

Glucose Tolerance Testing = TIGER TOP

There have been many glucose tolerance tests coming in without the tubes marked with the time or order of draw.

If we aren't provided this info the results will not be accurate.
You must label the tubes accordingly:

- Patients first and last name
- DOB
- Number the tubes 1-2 , 1-3 or 1-4

EXAMPLE:

1 hour glucose
#1 fasting
#2 1 hour draw

2 hour glucose
#1 fasting
#2 1 hour draw
#3 2 hours draw

3 hour glucose
#1 fasting
#2 1 hour draw
#3 2 hour draw
#4 3 hour draw